



Products & Capabilities

Alina[®]

User-friendly disposable pen injector platform
for variable and multi-dose treatments



Stevanato Group

MANAGING COMPLEXITY, DELIVERING VALUE



About Alina[®]

Alina[®] is an easy-to-use disposable pen injector platform with variable dosing.

The pen injector's performance, features and intuitive handling are compatible with established therapeutic regimens as well as innovator drug therapies that require a multidose delivery device such as diabetes and weight management.



Features

Suitable for **ISO 3mL glass cartridges** & optimized for Stevanato Group's range of glass solutions

Compatible with **standard** pen needles

Dose indicator, with registered design, for **improved readability**

Removable cap

Dial-back **dose error correction**



Benefits



**Ideally suited for
diabetes and weight
management**



**Visual and audible
feedback for dose setting,
correction & injection**




**Tested & verified
according to ISO11608
standards**




**Customization available
for primary packaging &
therapy needs**


Handling Steps




Remove cap




Attach the needle




**Hold the device
upright and prime**



**Turn the dial
to select the dose**



**Inject carefully
into a suitable
injection pad**




**Remove the needle
& recap the device**

Your Integrated Solutions Provider

Production of Alina® takes place at our state-of-the-art, ISO13485 certified, FDA-inspected manufacturing site in Germany.

We have invested in component tooling, automated sub-assembly equipment and optimized production processes to accelerate product launches and minimize up-front investments for customers.

Additionally, we have a range of optional products and services to further simplify your supply chain.





Headquarters

Via Molinella, 17
35017 Piombino Dese
Padova, Italy

**Global presence with operating units
and commercial offices in 9 countries.**

Visit our website to contact us.



stevanatogroup.com

ALINA_202402