

VIALS¹

| Vial Format | Nominal Capacity | Outer Diam. mm | Height mm. | Crimp Finish mm. | Secondary Packaging | No Vials x N&T/Tray | N&T/Tray x Kit | No of Vials x Kit | Design | Stopper dim. mm | Tot. Stopper x Kit | Alu cap dim. mm | Tot. Caps x Kit |
|------------------|------------------|----------------|------------|------------------|---------------------|---------------------|----------------|-------------------|--------|-----------------|--------------------|-----------------|-----------------|
| 2R ⁴ | 3 ml | 16 | 35 | 13 | Nest/tub | 100 | 6 | 600 | V9401 | 13 | 1000 | 13 | 1000 |
| | | 16 | 35 | 13 | Tray | 228 | 4 | 912 | V9401 | 13 | 1000 | 13 | 1000 |
| 6R ⁴ | 8 ml | 22 | 40 | 20 | Nest/tub | 48 | 6 | 288 | V9407 | 20 | 500 | 20 | 500 |
| | | 22 | 40 | 20 | Tray | 96 | 4 | 384 | V9407 | 20 | 500 | 20 | 500 |
| 10R ⁴ | 10 ml | 24 | 45 | 20 | Nest/tub | 48 | 6 | 288 | V9407 | 20 | 500 | 20 | 500 |
| | | 24 | 45 | 20 | Tray | 96 | 4 | 384 | V9407 | 20 | 500 | 20 | 500 |
| 20R ⁴ | 20 ml | 30 | 55 | 20 | Nest/tub | 24 | 6 | 144 | V9407 | 20 | 500 | 20 | 500 |
| | | 30 | 55 | 20 | Tray | 60 | 4 | 240 | V9407 | 20 | 500 | 20 | 500 |
| 25R ⁵ | 25 ml | 30 | 65 | 20 | Nest/tub | 24 | 6 | 144 | - | - | - | - | - |
| | | 30 | 65 | 20 | Tray | 60 | 4 | 240 | - | - | - | - | - |
| 30R ⁵ | 30 ml | 30 | 75 | 20 | Nest/tub | 24 | 6 | 144 | - | - | - | - | - |
| | | 30 | 75 | 20 | Tray | 60 | 4 | 240 | - | - | - | - | - |
| 50R ⁵ | 50 ml | 40 | 73 | 20 | Nest/tub | 16 | 6 | 96 | - | - | - | - | - |
| | | 40 | 73 | 20 | Tray | 28 | 4 | 112 | - | - | - | - | - |

¹ All vials are not blow back

⁴ Available with or without Add-on components

⁵ Available with Add-on-components upon request

RUBBER STOPPERS²

ALU CAPS³

² Compound: FM457

Coated: OmniFlex[®]

³ Disc Colour: Red

Type: Datwyler PrimeCaps[™]

SYRINGES

| Performance Level | Syringe Format | Inner Diameter | Outer Diameter | Cannula Length | Closure/ Needle Shield | No of Syringes x Nest/tub | N&T x Kit | No of Syringes x Kit |
|-------------------|----------------|----------------|----------------|----------------|------------------------|---------------------------|-----------|----------------------|
| Alba [®] | 0.5ml | 4.65 | 6.85 | - | ITC / OVS | 100 | 6 | 600 |
| | 1ml Long | 6.35 | 8.15 | 1/2" (12.7 mm) | RNS | 160 | 6 | 960 |
| | 2.25ml | 8.65 | 10.85 | 1/2" (12.7 mm) | RNS | 100 | 6 | 600 |
| Nexa [®] | 1ml Long | 6.35 | 8.15 | 1/2" (12.7 mm) | RNS | 160 | 6 | 960 |
| | 2.25ml | 8.65 | 10.85 | 1/2" (12.7 mm) | RNS | 100 | 6 | 600 |
| Fina [®] | 1ml Long | 6.35 | 8.15 | - | ITC | 100 | 6 | 600 |
| | 1ml | 8.65 | 10.85 | - | ITC | 100 | 6 | 600 |
| | 2.25ml | 8.65 | 10.85 | - | ITC | 100 | 6 | 600 |
| | 3ml | 8.65 | 10.85 | - | ITC | 100 | 6 | 600 |

ITC (Integrated Tip Cap) - RNS (Rigid Needle Shield)

CARTRIDGES

| Performance Level | Cartridge Format | Outer diameter mm | Height* mm | Finish mm | No Cartridges x N&T | N&T x Kit | No Cartridges x Kit | Configuration ** |
|-------------------|------------------|-------------------|------------|-----------|---------------------|-----------|---------------------|---------------------------|
| Nexa [®] | 1,5ml | 8.65 | 56.95 | 7.15 | 100 | 6 | 600 | Pre-capped with Lined cap |
| | 3ml | 11.6 | 62.3 | 7.15 | 100 | 6 | 600 | Pre-capped with Lined cap |

*The height refers to the glass cartridge without the assembled closure. **Pre-capped Cartridge (upside-down)